

# 96-Well Genomic DNA Kit (Plant)

*For research use only*

<b>Sample</b>	: up to 100 mg of fresh or 25 mg of dry plant tissue
<b>Yield</b>	: up to 30 µg/well
<b>Format</b>	: 96-well plates
<b>Operation time</b>	: within 60 minutes

**Geneaid**



[www.geneaid.com](http://www.geneaid.com)

## Introduction

The 96-Well Genomic DNA Kit (Plant) provides an efficient method for purifying total DNA (including genomic, mitochondrial and chloroplast DNA) from plant tissue and cells. Samples are disrupted by both grinding in liquid nitrogen and lysis buffer incubation. The lysate is treated with RNase A and then filtered to remove cell debris and salt precipitates. In the presence of the binding buffer coupled with chaotropic salt, the genomic DNA in the lysate binds to each column of the binding plate. The contaminants are removed with a Wash Buffer (containing ethanol) and the purified genomic DNA is eluted by a low salt Elution Buffer or TE. The procedure does not require DNA phenol extraction or alcohol precipitation and can be completed in 1 hour. The purified genomic DNA is ready for use in PCR, Real-time PCR, Southern Blotting and RFLP.

## Quality Control

The quality of the 96-Well Genomic DNA Kit (Plant) is tested on a lot-to-lot basis by isolating genomic DNA from a 50 mg young leaf sample. More than 10 µg of genomic DNA is quantified with a spectrophotometer and checked by electrophoresis.

### Kit Contents

Name	GPP02	GPP04	GPP10
GP1 Buffer	100 ml	200 ml	500 ml
GPX1 Buffer	100 ml	200 ml	500 ml
GP2 Buffer	25 ml	50 ml	125 ml
GP3 Buffer* (Add Isopropanol)	150 ml (300 ml)	150 ml (300 ml)	150 ml x 2 (300 ml x 2)
W1 Buffer	100 ml	200 ml	500 ml
Wash Buffer** (Add Ethanol)	50 ml (200 ml)	50 ml (200 ml)	50 ml x 3 (200 ml x 3)
Elution Buffer	30 ml	80 ml	100 ml
RNase A (10 mg/ml)	1 ml	2 ml	5 ml
DNA Binding Plate	2 pcs	4 pcs	10 pcs
0.35 ml Collection Plate	2 pcs	4 pcs	10 pcs

### Order Information

Product Name	Package Size	Cat. No.
Genomic DNA Mini Kit (Blood/Cultured Cell)	100/300 preps	GB100/300
Genomic DNA Midi Kit (Blood/Cultured Cell)	25 preps	GDI25
Genomic DNA Maxi Kit (Blood/Cultured Cell)	10/25 preps	GDM10/25
Genomic DNA Mini Kit (Tissue)	50/100/300 preps	GT050/100/300
Genomic DNA Mini Kit (Plant)	100 preps	GP100
Genomic DNA Maxi Kit (Plant)	10/25 preps	GPM10/25
96-Well Genomic DNA Kit	2/4/10 x 96 Wells	GBP02/04/10
96-Well Genomic DNA Kit (Plant)	2/4/10 x 96 Wells	GPP02/04/10
Vacuum Manifold (Accessories)	1 set	ZFV04

\* Add Isopropanol to the GP3 Buffer prior to initial use (see the bottle label for volume).

\*\* Add absolute ethanol to the Wash Buffer prior to initial use (see the bottle label for volume).

## Caution

The components contain irritants. During operation, always wear a lab coat, disposable gloves, and protective goggles.

## Note

Due to various plant species containing different metabolites such as polysaccharides, polyphenols, and proteins, we provide two different lysis buffers to satisfy the various plant samples. The standard protocol uses GP1 Buffer for lysis of plant samples. For most common plant species, the buffer system ensures purified DNA with high yields and high quality. Alternatively, GPX1 Buffer is provided with the kit. The different detergent in this lysis buffer is suitable for some plant samples with high polysaccharide content.

## 96-Well Genomic DNA Kit Plant Protocol

- Add Isopropanol to the GP3 Buffer prior to initial use (see the bottle label for volume).
- Add absolute ethanol to the Wash Buffer prior to initial use (see the bottle label for volume).
- Additional Requirements: centrifugation system for 96-well plates, 2 ml collection plates, Isopropanol, absolute ethanol

Step 1 Tissue Dissociation	<ul style="list-style-type: none"> <li>● Cut off 50 mg (up to 100 mg) of fresh or frozen plant tissue or 10 mg (up to 25 mg) of a dried sample.</li> <li>● Grind the sample under liquid nitrogen to a fine powder and transfer it to each well of a 2 ml collection plate (some plant samples can be disrupted without liquid nitrogen).</li> </ul>
Step 2 Lysis	<ul style="list-style-type: none"> <li>● Add <b>400 µl of GP1 Buffer</b> (or <b>GPX1 Buffer</b>) and <b>5 µl of RNase A</b> into each well of the 2 ml collection plate and mix by pipetting. <b>Do not mix GP1 Buffer (GPX1 Buffer) and RNase A before use.</b></li> <li>● Incubate at 65°C for 10 minutes. During incubation, gently shake the plate every 5 minutes. At this time, preheat the required <b>Elution Buffer</b> (200 µl per sample) to 65°C (for Step 5 DNA Elution).</li> <li>● Add <b>100 µl of GP2 Buffer</b> to each well of the 2 ml collection plate and mix by pipetting.</li> <li>● Incubate at 4°C for 3 minutes.</li> <li>● Centrifuge at 1,000 x g for 10 minutes.</li> <li>● Transfer the supernatant to a new 2 ml collection plate.</li> </ul>
Step 3 DNA Binding	<ul style="list-style-type: none"> <li>● Add a 1.5 volume of <b>GP3 Buffer</b> (Isopropanol added) to the lysate and gently shake immediately for 5 seconds (eg. add 750 µl GP3 Buffer to 500 µl lysate).</li> <li>● Place a <b>DNA Binding Plate</b> on a new 2 ml collection plate.</li> <li>● Transfer 700 µl of the mixture to each well of a <b>DNA Binding Plate</b>.</li> <li>● Centrifuge at 1,000 x g for 2 minutes.</li> <li>● Discard the flow-through from the 2 ml collection plate. Place the <b>DNA Binding Plate</b> back on the 2 ml collection plate and transfer the remaining mixture to each well of the <b>DNA Binding Plate</b>.</li> <li>● Centrifuge again at 1,000 x g for 2 minutes.</li> <li>● Discard the flow-through and place the <b>DNA Binding Plate</b> back on the 2 ml collection plate.</li> </ul>
Step 4 Wash	<ul style="list-style-type: none"> <li>● Add <b>400 µl of W1 Buffer</b> to each well of the <b>DNA Binding Plate</b>.</li> <li>● Centrifuge at 1,000 x g for 30 seconds.</li> <li>● Discard the flow-through and place the <b>DNA Binding Plate</b> back on the 2 ml collection plate.</li> <li>● Add <b>600 µl of Wash Buffer</b> (ethanol added) to each well of the <b>DNA Binding Plate</b>.</li> <li>● Centrifuge at 1,000 x g for 30 seconds.</li> <li>● Discard the flow-through and place the <b>DNA Binding Plate</b> back on the 2 ml collection plate.</li> <li>● Centrifuge again for 3 minutes at 1,000 x g to dry the column matrix.</li> </ul>
Step 5 DNA Elution	<p>Standard elution volume is 100 µl. If less sample is to be used, reduce the elution volume (30-50 µl) to increase DNA concentration. If higher DNA yield is required, repeat the DNA Elution Step to increase DNA recovery and the total elution volume to approximately 200 µl.</p> <ul style="list-style-type: none"> <li>● Transfer the dried <b>DNA Binding Plate</b> to a <b>0.35 ml Collection Plate</b>.</li> <li>● Add <b>100 µl of preheated Elution Buffer</b> or TE to the center of each well of the <b>DNA Binding Plate</b>.</li> <li>● Let stand for 3-5 minutes or until the <b>Elution Buffer</b> or TE is absorbed by the matrix.</li> <li>● Centrifuge at 1,000 x g for 1 minute to elute the purified DNA.</li> </ul>

## Troubleshooting

Problem	Possible Reasons/Solution
Clogged Column	<p><b>Too much sample was used.</b></p> <ul style="list-style-type: none"> <li>● Reduce the sample volume</li> </ul>
Low Yield	<p><b>Precipitate was formed at the DNA Binding Step</b></p> <ul style="list-style-type: none"> <li>● Reduce the sample material.</li> <li>● Prior to loading the plate, break up the precipitate in the ethanol-added lysate.</li> </ul> <p><b>Incorrect DNA Elution Step</b></p> <ul style="list-style-type: none"> <li>● Ensure that the Elution Buffer or TE is added to the center of each well of the DNA Binding Plate and is absorbed completely.</li> </ul>
Eluted DNA does not perform well in downstream applications	<p><b>Incomplete DNA Elution</b></p> <ul style="list-style-type: none"> <li>● Elute twice to increase yield.</li> </ul> <p><b>Residual ethanol contamination</b></p> <ul style="list-style-type: none"> <li>● Following the Wash Step, dry the DNA Binding Plate with additional centrifugation at 1,000 x g for 5 minutes or incubate at 60°C for 5 minutes.</li> </ul>